Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

COBRA Cooled™ Surgical Probe With Integrated Cable

Catalog Number 600-004*i*

U.S. Patent No. 5,651,780; 5,769,847; 5,797,905; 5,810,802; 5,871,523; 6,030,382; 6,106,522; 6,129,724; 6,387,092; 6,447,506;
6,471,699; 6,500,172; 6,142,994; 6,425,895; 6,120,496; 6,241,754; 6,416,505; 6,830,576; 7,115,122

Instructions for Use

TABLE OF CONTENTS

English	2
Dutch	6
French	10
German.....	14
Italian	18
Spanish	22

Sterile – Single Use Only

Estech, Inc.
2603 Camino Ramon Suite 100
San Ramon, CA 94583 USA
Phone: +1 925-866-7111
Fax: +1 925-866-7117
E-mail: info@estech.com

European Authorized Representative:
Medical Partners BV
Groenendijk 39
4926 RE Lage Zwaluwe
The Netherlands
Phone: +31.168.48.37.32
Fax: +31.168.48.47.90

ENGLISH

Instructions for Use

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these directions. Failure to do so may result in patient complications.

ESTECH relies on the physician to determine, assess and communicate to each patient all foreseeable risks of the procedure.

Description:

The ESTECH COBRA Cooled™ Surgical Probe with Integrated Cable (Probe) has a malleable distal section. These Probes are designed to conform to the specific anatomy of the tissue area to be coagulated.

The distal section of the Probe five 12.5mm electrodes spaced 2mm apart along the body. Any combination of coagulating electrodes may be used. The Probe is connected to the Electrosurgical Unit (ESU) by a cable permanently attached to the handle of the probe. Instructions for the use of the ESU and cable may be found in the ESU Operator's Manual.

For your convenience, a positioning tool is included in the packaging. This device is intended only to be used to help position and hold the probe against tissue.

Indications:

The ESTECH COBRA Cooled™ Surgical Probe with Integrated Cable is intended for the coagulation of cardiac tissue using radiofrequency (RF) energy during cardiac surgery. The Probe can be used during general surgery to coagulate soft tissues. The Probe may also be used to coagulate blood and soft tissue to produce hemostasis.

Contraindications:

- Local or systemic infection

Warnings:

- Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found call your ESTECH representative.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- The device(s) should be used by physicians thoroughly trained in the techniques of invasive surgical procedures.
- The operator should keep the set temperature and power limit as low as possible to achieve the desired end effect. This avoids excessive thermal damage to tissue, or collateral damage to adjacent tissue not intended for coagulation.
- Care should be taken to assure that the probe is not in contact with tissue other than that to be coagulated to avoid inadvertent tissue damage.
- Care should be taken when using the probe in proximity to vascular and nerve tissue to avoid inadvertent tissue damage.
- Care should be taken to thermally isolate the tissue to be coagulated when anatomically possible to avoid damage to unintended tissues or structures.
- Following RF coagulation, visual inspection of underlying tissues should be routinely performed to rule out the presence of inadvertent tissue damage.
- Care should be taken to assure that the probe is not in contact with other surgical instruments, staples or other objects while coagulating. Inadvertent contact with objects while coagulating could lead to conduction of RF Energy or heat and unintentional coagulation of tissues in contact with that object.
- Care should be taken when positioning surgical probe or positioning tool to prevent perforation or other damage to adjacent tissue.

**Precautions:**

- Before using, inspect for physical damage, including electrical insulation on the cable and the Probe shaft. Replace Probe if damaged.
- If using a TEE probe, care should be taken to withdraw the TEE probe prior to ablation to avoid compressing the esophagus against the left atrium during ablation.
- Many variables, including patient anatomy, pathology, and surgical techniques, may influence procedural outcomes. Patient and procedure selection is a responsibility of the medical professional.
- Avoid attaching the device or ablating over extremely thin, fragile, or aneurysmal tissue.
- Do not attempt to operate the system before thoroughly reading the ESU Operator's Manual and Instructions for Use. For future reference, keep these documents in a convenient, easily accessible place.
- While the distal portion of the Probe is designed to be malleable and/ or flexible, to conform to the anatomy of the area to be coagulated, excessive or rough shaping of the probe may damage internal components of the Probe.
- Pacemakers and implantable cardioverter/defibrillators can be adversely affected by RF signals. Refer to the manufacturer's Directions for Use.
- Dispersive Indifferent Patch (DIP) electrodes used with the Probe should be applied carefully according to the manufacturer's directions. Poor or incomplete contact of the DIP electrodes may result in skin burns. The use of DIP electrodes which meet or exceed ANSI/AAMI HF-18 requirements is recommended.
- The risk of igniting flammable gases or other materials is inherent in the application of RF power. Precautions must be taken to restrict flammable materials from the area where tissue coagulation is performed.
- Electromagnetic interference (EMI) produced by the ESU during the delivery of RF power may adversely affect the performance of other equipment.
- The Probe is intended for use only with the ESU and accessories.

Storage:

- Store in a cool dry place.

Inspection Prior to Use:

Carefully inspect the package prior to use for any breach of the sterile barrier or damage to the contents. If the sterile barrier integrity is compromised or the contents damaged, contact your ESTECH representative.

Directions for Use

1. Inspect the Probe, positioning tool and all packaging materials carefully. Open the package using sterile technique.
2. Ensure that DIP electrodes are in good contact, according to the manufacturer's directions.
3. There is a single cable which exits the Probe handle. The connector of the cable fits into the black receptacle on the face of the ESU. The black receptacle is located on the right hand side behind the sliding door. Gently shape the distal portion of the catheter to conform to the anatomy of the area to be coagulated. The positioning tool may be used to help position and hold the Probe against the tissue.
4. Attach one female luer fitting to a fluid pump set at 500ml/hr flow rate. Either sterile water or normal saline can be used as a coolant. The other female luer fitting should be attached to a receptacle to catch coolant as it exits the probe. The fluid pump should be operating whenever the ESU is supplying energy.
5. Select electrodes to be activated on the ESU.
6. Activate selected electrodes by depressing the yellow on/off switch on the ESU.
7. Radiofrequency energy may be discontinued by depressing the yellow on/off switch on the ESU.

Ablation Time (seconds) COBRA Cooled™ Surgical Probe at 40W/Electrode*

Tissue Thickness	70° C	80° C
5mm	60	60
10mm	120	90

*Recommended ESU Settings: 40 Watts/coil, 70°C, 60 seconds (endocardial);
40 Watts/coil, 70°C, 120 seconds (epicardial)

Complications

The following potential risks or discomforts may be associated with electrosurgical procedures. The frequency and severity of these events can vary, and may necessitate additional medical intervention, including surgery. Strict adherence to the forgoing instructions before use will help reduce the incidence of complications.

Allergic reaction, Arrhythmias, Cardiac or respiratory arrest, Cardiac valve damage, Chest pain, Damage to vessel intima or cardiac ultrastructures, Death, Hematoma / ecchymosis, Hemorrhage, Infarction, Infection, Perforation, Pericardial effusion, Pericarditis / pleuritis, Pseudoaneurysm, Pulmonary edema, Sinus or AV node injury, Stroke, Tamponade, Thrombosis, Vasovagal reaction.

Disposal

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Warranty and Limitations

ESTECH warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond ESTECH's control directly affect the instrument and the result obtained from its use. ESTECH's obligation under this warranty is limited to the repair or replacement of this instrument and ESTECH shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. ESTECH neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. ESTECH assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instrument.

For Technical Information Contact:

Endoscopic Technologies, Inc. dba Estech

2603 Camino Ramon Suite 100

San Ramon, CA 94583 USA

Phone: +1 925-866-7111

Fax: +1 925-866-7117

E-mail: info@estech.com

European Authorized Representative:

Medical Partners BV

Groenendijk 39

4926 RE Lage Zwaluwe

The Netherlands

Phone: +31.168.48.37.32

Fax: +31.168.48.47.90

This Product Contains No Latex.

This Product May Contain DEHP.

Graphic Symbols for this Device Labeling

<p>Catalogue Number</p> 	<p>Caution-See instructions for Use</p> 	<p>Manufacturer/Company Address</p> 	<p>Sterilized using Ethylene Oxide</p> 
<p>Batch Code</p> 	<p>Single Use Only</p> 	<p>Expiration Date</p> 	<p>Authorized Representative in the European Union</p> 
<p>Quantity</p> 	<p>Date of Manufacture</p> 	<p>Non-Pyrogenic</p> 	<p>Does not contain Natural Rubber Latex</p> 
<p>Contains no di (2-ethylhexyl) phthalate (DEHP) released from polyvinyl chloride (PVC)</p> 	<p>"Consult Instructions For Use"</p> 	<p>CE Mark with Notified Body - class II and III medical devices comply with EC Directive 93/42/EEC</p> 	<p>"Caution: U.S. Federal law restricts this device to sale by or on the order of a physician or other licensed practitioner."</p> <p>Rx ONLY</p>

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